

SWISSLOG CUSTOMER SUPPORT BULLETINS

TRANSLOGIC 4-INCH UNIVERSAL CARRIER & COMPONENTS

When used in conjunction with an additional containment procedure, the Universal Carrier provides a secure environment to send fragile items and other specialized materials.



Universal Carriers are constructed of clear polycarbonate with self-closing latches and a side-opening design to give the user full access to the carrier interior. Carriers are bi-directional and can be loaded into the station in either orientation.

4-inch Universal Carriers are available with a variety of options

including color and BacStop as well as frosted-end and half-black, which can be used in non-TransLogic systems.

BacStop, a silver ion antimicrobial agent, is embedded within the body of the carrier, offering effective protection against many classifications of bacteria, yeast and fungi.

Ordering Information for Carriers

Standard 4-inch Universal Carrier	075-1710-05
Clear with BacStop	075-1710-12
Frosted-End Carrier	075-1710-09
Frosted-End Black Carrier	075-1710-21
Frosted-End Carrier with BacStop	075-1710-15
Black Carrier	075-1710-08
Red Carrier with BacStop	075-1710-17
Red Carrier with BacStop and Frosted-Ends	075-1710-19
Yellow Carrier with BacStop	075-1710-18
Yellow Carrier with BacStop and Frosted-Ends	075-1710-20
Half-Black Carrier	075-1710-07
Half-Black Carrier with Frosted-Ends	075-1710-16
Half-Red, Half-Clear Carrier	075-1710-23
Half-Red, Half-Clear Carrier with Frosted-Ends	075-1710-25
Half-Yellow, Half-Clear Carrier	075-1710-22
Half-Yellow, Half-Clear Carrier with Frosted-Ends	075-1710-24

TransLogic 4-Inch Universal Carrier Components

Wear Band Replacement Kit



Hinge Pins



Latch Replacement Kit



Ordering Information for Carrier Components

Wear Band Replacement Kit
 Latch Replacement Kit
 Hinge Pins

(repairs 25 carriers) 076-1137-01
 (repairs 10 carriers) 076-1394-02
 (repairs 5 carriers) 076-1488-02

Due to ongoing improvement programs, Swisslog reserves the right of production or design change without notice or obligation.

For More Information in North America:

Swisslog Customer Support
 Email: orderentry@swisslog.com
 USA: (800) 396-9666 | Fax: (303) 373-7871
 Canada: (877) 294-2831 | Fax: (905) 629-2799

©2009 Swisslog
 CM# CS-401

WWW.SWISSLOG.COM

swisslog