



MEDPICKRX™
PICK-TO-LIGHT STORAGE & PICKING SYSTEM

HEALTHCARE
SOLUTIONS



Visit the MedPickRx web page.

swisslog

MedPickRx™ Storage and Picking System

MedPickRx is a high-density automated pick-to-light pharmacy system that optimizes storing and picking unit doses. MedPickRx is an inventory tracking and management program that supervises inventory locations within a custom-configured grouping of storage modules: workstations, drawers, carousels, and refrigerators.

MedPickRx provides all of the benefits of vertical carousels at a fraction of the cost.

MedPickRx Pick-to-Light Inventory Management System



Workstation
& Modular
Storage



High-Density
Drawer Module



Carousel



Refrigerator

Efficient, Compact and Safe

Patient orders or dispensing cabinet replenishment orders are downloaded into MedPickRx. Drugs are listed sequentially for picking from individually defined storage pockets. Each pocket has a unique identification allowing LEDs to quickly and efficiently direct an operator to its location. Barcode scanners verify the correct location and drug for picking.



Efficient Picking

Unlike aisles of static shelving, the clustered MedPickRx storage modules enable technicians to pick orders with fewer steps. Color-coded LEDs enable multiple techs to be picking simultaneously from the same modules. Refrigerators extend the efficiencies of pick-to-light technology to refrigerated products.




Inventory Management

MedPickRx picks oldest first, minimizing expiration loss. Par levels are based on historical usage patterns, overall and seasonal. MedPickRx eliminates multiple open bottles of the same item and the associated waste. MedPickRx provides extensive reporting and inventory analysis including an audit trail for every operator action.



High-Density Storage

Counter depth pull-out drawers maximize storage space in the pharmacy. Storage pocket locations are maintained by the inventory management software to ensure continued space optimization. High-density drawer modules utilize vertical space.




Patient Safety

Pick-to-light technology minimizes picking errors inherent with picking from static shelving. This allows the segregation of "sound-alike/look-alike" medications, eliminating alphabetized based storage. Barcode scan-in/scan-out verifies that the correct pocket and correct item is stored and picked. First DataBank imaging provides visual confirmation of the correct item.



For More Information in North America:

Swisslog Healthcare Solutions
Email: healthcare@swisslog.com
USA: 800.764.0300
Canada: 877.294.2831 | 905.629.2400

www.swisslog.com
 [SwisslogHealthcareSolutions](#)
 [SwisslogNA](#)
 [SwisslogHCS](#)

swisslog