



MEDROVER™
MOBILE DISPENSING CABINET

HEALTHCARE
SOLUTIONS



Watch the
MedRover™
Mobile
Dispensing
Cabinet Video

MedRover™

Finally...

Innovation like you've never seen before. It's a workstation on wheels, medication cart and dispensing cabinet in ONE compact unit.

Swisslog offers MedRover, a mobile drug cabinet which dramatically increases dispensing safety and time at the patient bedside, by reducing the amount of time spent acquiring the medications and supplies from various locations. MedRover's flexible storage options and lightweight, ergonomic design make it an ideal solution for hospitals looking to reduce or augment existing stationary cabinets. MedRover holds patient-specific medications, floor stock, and supplies nurses need in a convenient cart that can be brought to the patient bedside for administration and verification, significantly reducing the number of miles walked each shift, while increasing patient safety.

Features:

- > Bedside verification
- > Simple touch interface
- > Triple lock security
- > Biometric authentication
- > Locked supply drawer
- > Locking wall dock
- > Pull-out work surface
- > Up to 108 storage bins (multiple sizes)
- > Specific bins for patient meds and floor stock



Benefits for Nurses:

- > Reduces or eliminates frequent and disruptive trips to stationary medication cabinets
- > Small, light, maneuverable, and sized for the patient's bedside
- > Unique, electronically-controlled containers securely deliver only selected medications
- > Supports clinical/bedside software and barcode checking
- > Detailed auditing and reporting of medication administration events

Benefits for Pharmacists:

- > Reduces overall costs by tracking all medications to prevent loss
- > Semi-automated refilling of medication containers
- > Secure, cabinet-style dispensing of patient medications
- > Flexible computer platform accommodates bedside verification software
- > Accommodates various replenishment models

Software:

MedRover utilizes two unique software programs, RoverRX™ and RoverRN™.

RoverRX™ for pharmacy communication with the cart

- > Supervisory system for all picking and return operations
- > System management and reporting

RoverRN™ for cart operation

- > Cart-controlled 5-Rights
- > Scheduled/PRN medications at the bedside



PC and EMR Software

Wall Dock Station
> Security
> Battery charging

Biometric Authentication
> Log on witnessing

Pull-Out
Work Surface

Locked Supply Drawer

Up to 99 Secure
Storage Containers
> Triple lock security
> Individually locked
containers
> Controlled drug compliant
> Large items (PCA, IV, etc.)
> Floor stock containers
> Patient-specific containers
> Isolated respiratory care
containers

Ergonomically Designed
> Lightweight plastics
> No corners
> Easy to rotate in tight spaces





Swisslog - *Simply Better*

- > Visionary Product Innovation
- > Swift Return-on-Investment
- > Outstanding Product Performance
- > Superior Product Reliability
- > Excellent Customer Service and Support



For More Information in North America:

Swisslog Healthcare Solutions
Email: healthcare@swisslog.com
USA: 800.764.0300
Canada: 877.294.2831 | 905.629.2400

www.swisslog.com
 SwisslogHealthcareSolutions
 SwisslogNA
 SwisslogHCS

