

# SWISSLOG SPECIFICATIONS

# MEDROVER™ MOBILE DISPENSING CABINET



## Containers



### Small Container

36 Locations



**Dimensions:** 2.5 x 6.5 x 2.75 H  
Ideal for floor stock medications

### Medium Container

18 Locations - Transferrable



**Dimensions:** 5.25 x 6.75 x 2.75 H  
Ideal for patient-specific meds

### Large Container

9 Locations

IVs | PCAs | Bulk Items | Supplies




**Dimensions:** 4.75 x 7 x 12 H  
Ideal for large bulk items, supplies and IV Piggybacks

## Locking Wall Dock



Security | Charging

## Simple Touch Interface



## Specifications

Weight: 105 lbs  
Lithium polymer battery life: 9 hours  
Wheel base: 22 inches  
Storage containers: 63 of various size  
UL 60601-1 design

## Features:

- > Accommodates PC or Tablet
- > Night light
- > Auto log-off

### For More Information in North America:

Swisslog Healthcare Solutions  
Email: [healthcare@swisslog.com](mailto:healthcare@swisslog.com)  
USA: 800.764.0300 | Canada: 877.294.2831 or 905.629.2400

[SwisslogHealthcareSolutions](#) [SwisslogNA](#) [SwisslogHCS](#)

[WWW.SWISSLOG.COM](http://WWW.SWISSLOG.COM)