

## SWISSLOG SPECIFICATIONS

### TRANSLOGIC STATION BIN LINER FOR 6-INCH PNEUMATIC TUBE SYSTEM STATIONS



Swisslog's Station Bin Liner for the 6-inch TransLogic pneumatic tube system station provides durable, easy cleanup of the station bin interior.

### Features & Benefits

- > Smooth, non-porous surface provides quick cleanup of interior spills
- > Simple installation
- > Easily removable
- > Long-lasting seamless construction eliminates wear-prone carpet and padding
- > Fits all 6-inch TransLogic station housings manufactured after 1/2005. The bin can be used in older 6-inch stations with modification of the station housing. Contact Customer Support at (800) 396-9666 for details.



**Part Number:**  
**56679701**

The Station Bin Liner kit includes the station bin, station ramp cover, hardware and installation instructions.

*Due to ongoing improvement programs, Swisslog reserves the right of production or design change without notice or obligation.*

### For More Information in North America:

Swisslog Healthcare Solutions  
Email: [healthcare@swisslog.com](mailto:healthcare@swisslog.com)  
USA: (800) 764-0300  
Canada: (877) 294-2831 | (905) 629-2400

©2009 Swisslog  
CM# PTS-211

**WWW.SWISSLOG.COM**

**swisslog**