

## SWISSLOG SPECIAL APPLICATIONS

## ALERT MESSAGING FOR TRANSLOGIC PNEUMATIC TUBE SYSTEMS



Swisslog's Alert Messaging provides communication to users and/or operations on key pneumatic tube system transport events—increasing user productivity, improving system maintenance, and reducing the potential for system down-time.

### Automated System Notifications

With Swisslog's TL 2007 pneumatic tube system control software and more recent software releases, hospitals can benefit from an optional Alert Messaging package, which offers important and immediate system event notification.

Alert Messaging provides both hospital department users and the hospital engineering staff with a dynamic communication tool monitoring PTS transport events. Notification using several facility communication products is available to report system interruptions or failures in the following categories:

- > Communications
- > User-related situations
- > Maintenance issues
- > Equipment failures

The Alert Messaging feature is compatible with the following facility communication products:

- > Pagers
- > PDA
- > Mobile phones
- > E-mail
- > Walkie-talkie phones (Nextel®)
- > And more

\*Nextel is a registered trademark of Nextel Communications and has no affiliation with Swisslog.

## Alert Messaging for TransLogic Pneumatic Tube Systems

Notification categories can be associated with various hospital departments. A typical "user event" could be the notification of a full station bin, preventing the arrival of additional carriers to that station.

An example of a "maintenance event" could be notification to engineering of a system component requiring a response, such as a transfer unit or blower package. This email would provide identification of the system component and a brief text descriptive message.



Obtaining Swisslog's Alert Messaging package provides more advanced communication features, such as wireless, voice, email and pager communication as well as system maintenance and scheduling control.

*Due to ongoing improvement programs, Swisslog reserves the right of production or design change without notice or obligation.*

### For More Information in North America:

Swisslog Healthcare Solutions  
Email: [healthcare@swisslog.com](mailto:healthcare@swisslog.com)  
USA: (800) 764-0300  
Canada: (877) 294-2831 | (905) 629-2400

©2009 Swisslog  
CM# PTS-501

[WWW.SWISSLOG.COM](http://WWW.SWISSLOG.COM)