



SWISSLOG HEALTHCARE

Automated Packaging and Dispensing  
for Correctional Facilities

# Transforming Medication Management for Corrections

Pharmacy staff at Correctional Facilities face numerous challenges when managing, preparing and administering medications. Balancing the legal rights of inmates while adhering to regulatory requirements and mitigating the impact of rising drug costs is driving these facilities to seek alternatives to traditional medication management models.

The InSite System transforms the process of medication dispensing in corrections facilities from a tedious, error prone set of tasks to an efficient, automated workflow. This ensures adherence to strict regulatory requirements as well as the legal rights of inmates. The InSite "on-demand" dispensing model is well positioned to meet the unique challenges of transient inmate populations.



“ Previously it took 4 nurses, in this facility alone, 3 hours each to set up their medications and then to go pass their medications. Now, with the medication packager, we can deliver it for my entire facility in less than 30 minutes.”

- San Bernardino County Sheriff's Department

## Advantages

- **Accuracy** – Eliminates unlabeled doses to minimize the risk of medication errors
- **Efficiency** – Dramatically reduces the time to prepare for a med pass by packaging up to 60 doses per minute in configurable, patient-specific pouches
- **Responsiveness** – Provides real time location information through a direct interface to the inmate tracking system, so staff can easily suspend a patient's medication while they are off-site, in transit or at court
- **Convenience** – Easily dispenses "keep on person" or "self-carry" medications to keep nursing and pharmacy focused on managing high value drugs and caring for inmates
- **Security** – Restricts access to stored medications with password-protection and locking metal doors


## How On Demand Packaging Minimizes Medication Waste

Wasted medications cost corrections facilities hundreds of thousands of dollars every year. Due to the transient nature of incarceration, such as inmates being transferred to other facilities, paroled or released, many prescriptions go unused. With traditional 30-day blister cards, the inmate's remaining medications must be destroyed to ensure regulatory compliance.

Today, forward-thinking facilities are taking back control of their pharmacy inventory by packaging and dispensing medications on demand. The automated system tracks a patient's current location based on the latest data – even if they are off-site or at court. That means only the medications needed for a specific med pass are dispensed.

# InSite® In-Facility Medication Packaging and Dispensing System

The InSite System stores up to 240 different medication types in a single medication dispensing unit, providing automated, on-demand compliance packaging for oral solid medications prior to administration.



## Medication Security

Restricts access to stored medications with password protection and locking metal doors

## High-Speed Packaging

Packages up to 60 doses per minute with a broad variety of sorting options

## Access & Efficiency

The average InSite canister holds about 300 doses and can be filled and processed in roughly the same time it takes to create a single 30-day blister pack

## Safety

Supports an Automatic Canister Recognition System (ACRS) for added convenience and packaging accuracy

## Inventory Management

Allows automatic inventory monitoring as well as lot and expiration tracking

## Training Videos

Concise "how to" videos help staff use the machine independently


## How It Works

With InSite, medication orders are transmitted to the pharmacist for clinical review. After approval, the medications are immediately available at the corrections facility and dispensed from the InSite System with multi-dose, patient-specific packaging.

- Before delivery, each canister is filled, security-sealed, verified and signed by the pharmacist
- Medications are dispensed only as needed and for a specific patient.
- Changes in prescribed medications or patient residency will no longer produce unused blister packs that need to be destroyed or returned.
- InSite provides inventory and dispensing information for down-to-the-dose accountability at each facility



The InSite System is deployed to over 25,000 patients in corrections centers across the country, packaging over 32 million doses annually.



## Partnering with Swisslog Healthcare

Talyst, a market leader in pharmacy automation solutions, is now part of Swisslog Healthcare. Together, we're transforming medication management for long-term care facilities and pharmacies.

At Swisslog, we strive to lead change for better care. At the core of this vision is a focus on improving workflows and reducing the time clinicians spend doing routine and repetitive tasks—enabling more time to care for patients and residents. In addition to long-term care automation, our solutions and services extend across the continuum of care, including transport, medication and supply chain automation for hospitals and health systems.

Headquartered in Buchs, Switzerland, Swisslog Healthcare has been delivering best-in-class automation solutions for more than 100 years. We are a global organization with 2,500 team members in more than 20 countries and clients in more than 50 countries.



*InSite is a registered trademark of Swisslog AG.  
Swisslog systems may be covered by one or more patents. See [swisslog.com/patents](http://swisslog.com/patents) for details.*

Swisslog Healthcare  
healthcare.us@swisslog.com  
800.764.0300  
Canada: 877.294.2831 | 905.629.2400

[swisslog.com/corrections](http://swisslog.com/corrections)

**swisslog**  
healthcare